Full-Leg WoundRAP™ Instructions

- 1. Bring both zipper pulls to the top of the leg wrap.
- 2. Hold the left side of the zipper and pull down on the right side to separate the zipper, opening the leg wrap.
- Holding the wrap from the top, place the wrap around the leg, letting gravity naturally help the length of the wrap hang down.
- 4. Pull the top of the wrap around the leg, bring the ends together and connect the zipper.
- 5. Close the zipper to the bottom of the wrap. You will have to pull the wrap together as you gradually close the zipper. It is a compression bandage and it may take some effort to get the zipper closed at first, especially if there is swelling in the leg, but the soft compression is actually good for swollen wounds.
- 6. Secure the top elastic band by bringing it around the leg and through the D-ring, then Velcro it to itself. Give enough tension to secure the wrap but be careful not to over tighten, as this will cause it to become a tourniquet cutting off circulation. Wrap bottom elastic band around the pastern and secure with Velcro.
- 7. The two—way zipper makes it easier to check a wound and change out the dressings—from the top down or the bottom up—without removing the wrap completely.















