


## Foot Support Foam Pad

### TO REMOVE/REPLACE:

1. Remove old pad from **Foot Support**.
2. Clean the surface of the support.
3. Place the curved pad on the  
Foot Support and slide on (Figure 1).



*For more information or to re-order,  
you can visit [www.sunmedica.com](http://www.sunmedica.com) or give us a call  
and we will be happy to assist you!*



[www.sunmedica.com](http://www.sunmedica.com)

P 530.229.1600 | 1.800.995.8715 | F 530.229.9457

1661 Zachi Way, Redding, CA 96003